




Flex Coil Accessories

The optional Stabilizer and Pad Kits add extra support and comfort when additional stability or immobilization is required for the foot/ankle, knee or elbow. An integrated mechanism allows attachment of the coil to the base, eliminating additional pads or straps. The Flex Coil Pediatric Stabilizer enables high resolution imaging by combining two 1.5T GEM Flex Antenna arrays.




Stabilizer Foot/Ankle

Provide stable and comfortable support with a simple, lightweight, and adjustable stabilizer to capture foot or ankle imaging.




Stabilizer Elbow/Wrist

Comfortable support that is simple to use for both elbow and wrist imaging.



Stabilizer Knee

A simple to use stabilizer that avoids using hard plastics to provide additional knee imaging comfort.




Pediatric Stabilizer

Maximize your investment by adding 32-channel pediatric imaging. Each coil features 16 elements optimized to provide deep penetration SNR that is as good as, or better than, comparable rigid coils. The Pediatric Stabilizer supports the 1.5T GEM Flex Large coil to wrap around the head providing form-fitting comfort and immobilization. The 1.5T GEM Flex Medium lays flat on Pediatric Stabilizer for spine imaging. Depending on the orientation, parallel imaging acceleration factors of 2 or higher are practical.

Pad Kit

Add additional comfort where patients need it most with a unique set of pads to get the most comfortable fit while using the stabilizers.



Accessory Specifications

Component	Dimensions	Weight
Pediatric Stabilizer Base	70.1 cm x 41.4 cm x 16.8 cm	3.2 kg
Stabilizer Foot/Ankle without Slide	34.0 cm x 34.0 cm x 25.7 cm	1.8 kg
Stabilizer Knee	35.6 cm x 25.4 cm x 25.4 cm	0.7 kg
Stabilizer Elbow/Wrist	24.9 cm x 16.3 cm x 15.7 cm	0.5 kg
Pad Kit	Sizes Vary	Sizes Vary

Ordering Information

Item	Catalog
Pediatric Stabilizer	NC074003-D
Stabilizer Foot/Ankle Kit	NC074000-D
Stabilizer Knee Kit	NC074001-D
Extremities Positioner Kit	NC074004-D
3.0T Pad Kit	NC026009

For more Information

Email: sales@neocoil.com
Telephone: 888-MRI-COIL (674-2645)
www.neocoil.com